


Lab instrument log acquisition and analytics

Build a data pipeline to automate lab instrument log ingestion, collate metadata, and transform and downstream to FDA compliant electronic laboratory notebooks. Eliminate multiple second-person reviews for a [drug discovery and development](#) process. Experience easy search capabilities and REST API support.



- All instruments create the output log file in their respective network file system (NFS) shares on the machine attached to the instrument.
- An agent (virtual machine) is deployed in the user environment, and is used to read data from their NFS storage.
- An **AWS DataSync** task created in the customer's AWS account syncs data from NFS to a customer-owned **Amazon S3** bucket.
- A lab scientist invokes a client application to look for an experiment output log file in **Amazon S3** for an associated experiment ID.
- The client application gets a token from **Amazon Cognito** and invokes the API with the experiment ID and the **Amazon S3** file path (and possibly more IDs as needed from the client application).
- Amazon API Gateway** validates the token from **Amazon Cognito**, then calls the processing **AWS Lambda** function (one **Lambda** function for each type of log file), and returns response after the **Lambda** function runs.
- The **Lambda** function reads the file from **Amazon S3**, and adds the experiment ID as a tag to the **Amazon S3** log file.
- The **Lambda** function transforms the **S3** file, applies any business rules, converts the log to one or multiple JSON records, and pushes data to **Amazon Elasticsearch Service**.
- If the processing **Lambda** function fails, it writes the **S3** file info to an error monitoring queue. A **Lambda** function based on a pre-defined timed event can read from the queue to reprocess or to inform support staff.
- Optionally, a processing **Lambda** function maintains the **S3** file info with metadata (experiment IDs) in an **Amazon DynamoDB** table which can be consumed by other applications.
- The processing **Lambda** function sends the JSON response to **Amazon API Gateway**, which is posted to an electronic lab notebook (ELN) via Mobile Client.
- The lab scientist can access data in ELN, and can also discover and query data via Kibana for other experiments.

